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Sample preparation guidelines for MALDI-TOF

In MALDI-TOF-MS analysis samples are first dissolved and dispersed in a matrix (usually an organic acid), placed on a MALDI plate and, after evaporation of the solvent, the sample is introduced into the MALDI source where the sample-matrix crystals are irradiated by a laser beam.

The concentration of the sample to be delivered depends on the type of analyte, taking already into account the dilutions required with the matrix. Thus:

1 - For polymers, organic compounds and complexes the amount of sample needed is about 1 mg. The sample is dissolved in the solvent indicated by the applicant, which is preferably miscible with the solvent matrix.

2 - For peptides and proteins in solution, the concentration depends on molecular weight. For proteins with a mass below 10 kDa are needed about 10 to 20 pmol/ μ l. For proteins with masses exceeding 10 kDa the quantity to be delivered must be 100 to 200 pmol/ml. The proteins (or a mixture of peptides) may also be delivered lyophilized.

3 - Samples of protein from electrophoresis separated on polyacrylamide and stained with *Coomassie* blue or silver (silver compatible with mass spectrometry – please, contact the technician), can be delivered on the same gel hydrated. In the case of separate bands, the band/protein stain should be trimmed to the minimum possible gel and placed in eppendorf tube covered with type I water. It is essential to send a copy or a scanner image of the gel along with samples. It is strongly recommended the use of powder-free gloves for sample manipulation for identification by PMF, in order to avoid contamination with keratins.

Reagents compatible with MALDI-TOF/MS:

Always use reagents with HPLC purity and type I water.

- A) The compounds that do not interfere with the MALDI-TOF are: trifluoroacetic acid, formic acid, beta-mercaptoethanol, DTT, volatile organic solvents, chloridric acid, ammonium hydroxide, acetic acid.
- B) Compounds with concentrations below 50 mM do not interfere with mass spectrometry: HEPES, MOPS, Tris, NH₄OAc, octilglucoside.
- C) C) To avoid the following compounds: glycerol, sodium azide, DMSO, SDS, phosphates, sodium chloride, 2M urea, 2M guanidine.